


INSTRUCTION FOR HISTOPATHOLOGY SAMPLING



Objective

Describe the correct sampling of tissues for processing and subsequent histological diagnosis.

Materials and supplies: Scalpel handle, Scalpel blades, Scissors, Anatomic and surgical tweezers, Fixing solution: 10% buffered formalin, Container and Indelible marker pen.

Development of the process

- Select the culture units (cages or tanks) and prepare a smooth surface for sampling (Image 1).
- Label the containers with the fixing solution (10% buffered formalin) (Image 2). 10 parts of solution per 1 part of tissue are required (Image 3).
- Samples should be obtained only from fresh tissue. Those samples that show the minimum degree of autolysis **DO NOT HAVE ANY DIAGNOSTIC VALUE**.
- Scalpel cuts must be "clean", avoiding tearing the tissue. It is not recommended to use a rat tooth tweezers (Image 4).
- Sample size:** 0.5 cm thick (Image 5).
- Fish less than 4 cm in length, introduce them completely in the fixative, making an incision along the abdominal cavity without damaging the digestive tract.
- Pyloric caeca:** make a longitudinal cut and take the anterior portion to obtain pancreatic tissue.
- Kidney:** Take preferably medium kidney.
- Gills:** lift the operculum, select first branchial arch and cut including a piece no larger than 0.5 cm. **Injury:** the sample must include normal and damaged tissue.
- Eyes:** the sample must include the complete eyeball.
- Keep in a cool place between 1-20 ° C. Use gel pack if necessary. Do not expose to direct sunlight. Do not freeze.



Image 1: Materials for sampling



Image 2: Labeling of samples




Image 3: Fixing solution (10% buffered formalin, 1:10 (tissue:solution))




Image 4: Clean tissue cutting




Image 5: Sample size, 0.5 cm thick